


# VitaSIRO solo<sup>TM</sup>

## CT/NG/TV/MG Assay




### Detect and Defend with Cutting-Edge Technology

In a world where pathogens lurk undetected, safeguarding your health becomes paramount. Meet our revolutionary solution – VitaSIRO solo<sup>TM</sup> CT/NG/TV/MG Assay.

Our product utilizes state-of-the-art microfluidic channels to extract nucleic acids swiftly and efficiently, coupled with RT-qPCR technology for precise identification of Chlamydia Trachomatis (CT), Neisseria Gonorrhoeae (NG), Trichomonas Vaginalis (TV), and Mycoplasma Genitalium (MG) from urine or vaginal swab samples.



### Advantages:

- 
**Quick & Convenient Usage**  
 Ready-to-use cartridge with only 1 minute of hands-on time.
- 
**Sensitive & Rapid Results**  
 Sample to RT-qPCR results within 50 minutes, enabling testing and treatment within 1 hour.
- 
**Easy Storage & On-Demand Testing**  
 Room-temperature storage for easy access and on-demand testing to lessen anxiety.

### Ordering Information

Product Name	Product Code	Specification
VitaSIRO solo <sup>TM</sup> Instrument	PCRAC1101	1
VitaSIRO solo <sup>TM</sup> CT/NG/TV/MG Assay	PCRAD1107 PCRAD1117	12 Test / box 12 Test + 12 Dropper / Box

#### Credo Diagnostics Biomedical Pte. Ltd.

491B River Valley Road #06-01 Valley Point Singapore 248373

✉ [service@credodxbiomed.com](mailto:service@credodxbiomed.com)

🌐 [credodxbiomed.com](http://credodxbiomed.com)

🌐 [linkedin.com/company/credo-diagnostics](https://www.linkedin.com/company/credo-diagnostics)