

# VitaSIRO solo<sup>TM</sup>

## Human Papillomavirus (HPV) Assay

In a world where hidden threats to health persist, early detection is crucial. Meet our revolutionary solution - VitaSIROsolo<sup>TM</sup> HPV Assay. Our innovative solution employs advanced microfluidic channels for rapid nucleic acid extraction, combined with real-time PCR technology for accurate identification of high-risk HPV strains.

With a simple sample-in, result-out process, the HPV Assay is designed to work seamlessly with cervical specimens collected in ThinPrep<sup>®</sup> PreservCyt<sup>®</sup> Solution, ensuring precise and reliable results. This empowers healthcare providers to make informed decisions and enhance patient care.

### Advantages:



#### Quick & Rapid Results

Ready-to-use cartridge with only 1 minute of hands-on time and results reported in 45 minutes.



#### Inclusivity

Targets 14 high-risk HPV types with genotyping of HPV 16 and 18/45.



#### Portable & Adaptable

Small footprint and lightweight, perfectly fitting into various scenarios.

#### Target of Detection

- HPV 16
- HPV 18/45
- Other high-risk types:  
HPV 31, 33, 35, 39, 51, 52, 56, 58, 59, 66, 68

#### Specimen Type

- Cervical specimens collected in ThinPrep<sup>®</sup> PreservCyt<sup>®</sup> Solution



### Ordering Information

Product Name	Product Code	Specification
VitaSIRO solo <sup>TM</sup> Instrument	PCRAC1101	1 Unit / Box
VitaSIRO solo <sup>TM</sup> HPV Assay	PCRAD1104 PCRAD1114	12 Tests/ Box 12 Tests + 12 Droppers / Box

#### Credo Diagnostics Biomedical Pte. Ltd.

491B River Valley Road #06-01 Valley Point Singapore 248373

✉ [service@credodxbiomed.com](mailto:service@credodxbiomed.com)

🌐 [credodxbiomed.com](http://credodxbiomed.com)

🌐 [linkedin.com/company/credo-diagnostics](https://www.linkedin.com/company/credo-diagnostics)